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## TRACHEOSTOMY TUBES - INSTRUCTIONS FOR USE

These instructions for use refer to the following products:

REF: 31-XXXX, 32-XXXX, 33-XXXX, 34-XXXX, 35-XXXX, 36-XXXX, 37-XXXX, 38-XXXX, 39-XXXX

### Pictograms

|  |                                                               |  |                              |  |                                   |
|--|---------------------------------------------------------------|--|------------------------------|--|-----------------------------------|
|  | Catalogue number                                              |  | Manufacturer                 |  | Do not resterilize                |
|  | Batch code                                                    |  | Latex-free                   |  | Free from phthalates (incl. DEHP) |
|  | Single sterile barrier system<br>Sterilized by ethylene oxide |  | Consult instructions for use |  | CE marking with Notified Body no. |
|  | Use-by date                                                   |  | Medical device               |  | Date of manufacture               |
|  | Keep dry                                                      |  | Keep away from sunlight      |  | Store at 5-40 °C                  |
|  | Do not use if packaging is damaged                            |  | Do not re-use                |  |                                   |

Tracheostomy tubes are intended only for one patient. The maximum usage time of one tube should not exceed 29 days. It is recommended that the tubes be replaced frequently, as recommended by the doctor.

### Description

Assortment of sterile, single-use tracheostomy tubes, are characterized by the following properties:

- A low-pressure sealing cuff.
- A pilot balloon with a one-way valve.
- 15 mm connector
- The flange of the possibility of repositioning.
- Obturator for easy insertion.
- They are sterile until the package is opened or damaged.

1. The tube shaft is manufactured of transparent PVC. Each tube has a radiation line for confirmation of correct position of the tube using image diagnostics.
2. The cuff easily adopts to the natural shape of trachea and stops up the tracheal lumen also in the event of low pressure inside of the cuff. It is characterized by high residual volume (i.e. the quantity of air required to make impossible unrolling without causing increased pressure inside of each cuff). In such a way the trachea is sealed even upon low pressure in the cuff without bringing undesirable pressure to the mucosa of trachea.
3. If the device is used in a correct manner and according to the instructions, damage of mucous membrane or cartilage may be considerably limited. The pilot balloon indicates correct inflation of the cuff, though it does not reflect actual pressure in the cuff.
4. The olive end of obturator closes distant end of the tube securing in such a way non-traumatic placing of the tube in the tracheostomy opening.
5. The cuff supporting of the tube shaft, prevents the tube from dislocation.
6. The tube may possess 15 mm standard connector.
7. A cotton strap provides optimal fastening of the cuff on the patient's neck.

### Indications

The tracheostomy tubes are intended for patients after tracheotomy.

Stabilization of the tracheostomy after laryngectomy or tracheotomy; keeping the tracheostomy open.

In case of airway obstruction, after facial injuries with significant swelling, after injuries to the neck, chest or laryngeal burns. When using prolonged mechanical ventilation. During frequent suctioning of a significant amount of secretions from the trachea and bronchial tree. For protection against the aspiration of foreign bodies into the trachea, usually in unconscious patients or with impaired laryngeal reflexes. For patients who require continuous use of the tube.

### Contraindications

No absolute contraindications. Relative contraindications: allergy to components of the medical device, skin infection at the insertion site, coagulation disorders

### Instructions for use

1. Prior to set usage, it is essential to confirm its completeness and proper functionality. Verify the tubes with cuffs for any damages and ensure the air-tightness of the inflation system.
2. The procedure of tracheostomy should be performed by a medical professional.
3. One should select a proper size of the tube being guided by adopted medical procedures.
4. Gently remove the sterile tracheostomy tube from the package.
5. Prior to use check the cuff and testing balloon as well as valve of each tube inflating the cuff with air. To place the end-part of Luer syringe in housing of the valve and give an injection of air in the quantity sufficient to completely fill the cuff.
6. To check if the system is air-tight, check the pressure drop after 5 minutes
7. In order to prevent over-stretching of the cuff's membrane, the inflation pressure of the cuff during the test of air-tightness should not exceed 100 mbar.
8. After inflation test of the cuff, empty it completely of air.
9. Before inserting of the tube into the tracheostomy opening, one should make sure whether obturator is placed in such a way, that its tip closes the distal end of the tube. Keep the device in proper position using your thumb during the whole inserting process of the tube. Intubation may be facilitated by prior moistening of the ends of tube and obturator with a bit of lubricating water-soluble agent.
10. After inserting the tube into the trachea and checking correct positioning, remove the obturator.
11. Make sure, that the tube shaft is sitting in the tracheal opening without excessive force.
12. In case of a tracheostomy tube with adjustable flange, determine appropriate length of tube by moving the flange and lock it in place by twisting the cap
13. Aseptic gauze dressing should be put under the tube's cuff.
14. Tie the straps attached to the flange tightly enough so that the tube does not fall and not to exert pressure on external cervical veins
15. You should avoid disinfecting agents containing alcohol while localizing the tracheostomy opening.
16. You should never moisten gauze swabs with disinfecting agents containing alcohol.
17. You should inflate the cuff with a minimal quantity of air sufficient to ensure efficient tightness.
18. Hissing noises in the cuff area, especially when exhaling, indicate that the cuff does not completely fill the trachea, which is therefore not sealed. Ideally, a pressure of 15 mbar should be sufficient.
19. Remove the syringe from the housing of inflation valve after inflation of the cuff. Leaving of the syringe will cause opening of the valve and deflation of the cuff.
20. After intubation mount tube in place using tube holder, to prevent tube movements and accidental extubation.
21. Pressure in the cuff usually changes somewhat after a certain time due to diffusion of gases through the walls of cuff. Hence, we advise regular monitoring of the cuff's pressure.
22. Before deflation of the cuff you should remove by sucking secretion from the mouth, throat and trachea above the cuff.
23. Throw away the tracheostomy tube.

### Precautions

1. A patient suffering from tracheostomy can use the tracheostomy tube only if that decision has been taken by a physician.
2. The attending physician shall select an adequate type and size of the tube.
3. Insertion of the tracheostomy tube into the trachea always increases air resistance.
4. If there appears breathlessness, the patient should alert a nurse or a physician as soon as possible. In most cases it is efficient expectoration of secretion or mechanical removal by suction of air ways.
5. Under no circumstances should the end of the tube be closed with a phonation valve or decannulation plug or similar object while the cuff is inflated.
6. Closing of the tube (upon deflated cuff) using a decannulation cap, in order to prepare the patient for removal of the tube, may be performed only under supervision and at responsibility of a physician.
7. The cuff should be protected against contact with sharp edges.
8. One should avoid contact with coagulation electrodes or a laser beam, because PVC will emit toxic smoke, or it will cause ignition of surrounding environment enriched with oxygen.
9. The device is not proper for patients during radiotherapy.
10. Only SUMRI tubes can be used in MRI or CT environments.
11. Under no circumstances should the components of the set continue to be used if the shaft of the tube or the cuff is damaged (e.g. cracks on the surface, leaking cuff) and if the joints between the flange and the shaft are loosened, or the flange separates from the shaft
12. Discontinue the use of the tubes as soon as you notice signs of allergy

### Safety during imaging examinations (MRI/CT)

All tracheostomy tubes with cuff manufactured by SUMI, except for SUMRI models, are equipped with a valve containing a metal spring. Such tubes must not be used during magnetic resonance imaging (MRI) or computed tomography (CT) examinations. The presence of a metal spring in the valve or metal reinforcement in the tube wall may pose health risks to the patient and lead to inaccurate diagnostic results. If a tracheostomy tube is required during MRI or CT examination, it is recommended—in consultation with the attending physician—to use the SUMRI model, which does not contain any metal components.

### ATTENTION!

SUMRI tracheostomy tubes are safe for use during magnetic resonance imaging (MRI) and computed tomography (CT).

### Warnings

Current accepted medical procedures and techniques should be adhered to due to the risk of respiratory tract infection mainly with drug-resistant nosocomial strains of bacteria. Resterilization or reuse may pose a risk to the patient. The tube is for single-use only.

### Complications

Adhering to current approved medical techniques and having the procedure performed by qualified medical personnel helps reduce the risk of complications such as: tracheal and laryngeal stenosis, wound bleeding or infection, pneumothorax, mediastinal emphysema, tracheoesophageal fistula, blocking of the tracheostomy tube, inflammation and decubitus ulcer in the area of trachea, air leakage around the tracheostomy tube, respiratory tract infection, breathing difficulties and breathing effort, muscle lesions, laryngeal hemorrhage.

### Notes

1. The tubes should be stored in room conditions.
2. Do not use in case of damage of the package.
3. Do not use after the use-by date.

### Storage

Sterile tubes should be stored in a dry environment, protected from dust, sunlight, at a temperature of 5-40°C, in a way that protects them from moisture, direct sunlight, acid and organic solvent vapors and other harmful factors. During transport it is necessary to provide protection against weather conditions.

### Disposal

After use, dispose of this device in accordance with applicable national regulations on medical waste disposal.

### Responsibility – legal guidelines

The tracheostomy tubes manufactured by SUMI are subjected to physical and chemical testing and bio-compatibility testing. This allows their safe use and ensures highest quality. Carefully read the instructions for use before using the medical device. Tracheostomy tubes made by SUMI are intended for use only by qualified medical personnel and trained patients (users). The SUMI company will not be liable for the consequences caused by improper selection of tube size or by misuse of the tubes contrary to its intended purpose or provided instructions for use.

In the event of a complaint or an occurrence of a medical incident, keep the medical device with its primary packaging. It is also recommended to create photographic documentation. Every Medical Incident must be reported to the Manufacturer and relevant National Competent Authority.

SUMI sp. z o.o. sp. k.  
ul. Drobniarska 35, 05-070 Sulejówek, POLAND  
tel./fax. (+48 22) 783 30 23. www.sumi.com.pl

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