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Ed. 01 15.05.2024
IFU/RD/03

BRONCHIAL TUBE – INSTRUCTIONS FOR USE

Pictograms

	Catalogue number		Manufacturer		Free from phthalates (incl. DEHP)
	Batch code		Unique device identifier		Medical device
	Date of manufacture		Free from latex		Do not re-use
	Use-by date		Consult instructions for use		Keep dry
	CE marking with Notified Body no.		Keep away from sunlight		Temperature limit 5-40°C
	Sterilized using ethylene oxide, Single sterile barrier system				Do not use if packaging is damaged

These instructions for use refer to the following products:

21-xx10	21-xx60	21-xx61	22-xx10	22-xx60	22-xx61
21-xx20	21-xx70	21-xx71	22-xx20	22-xx70	22-xx71
21-xx22	21-xx80	21-xx81	22-xx22	22-xx80	22-xx81
21-xx41, 21-xx45	21-xx90	21-xx91	22-xx41, 22-xx45	22-xx90	22-xx91

- If an bronchial tube was lubricated before use, it one should check whether the lubricant has got into and clogged up lumen of the tube making in such a way ventilation impossible.
- One should periodically check position of the tube (auscultation, endoscopy of the bronchia, radiography).

Warnings

Because of a risk of infection of airways with some hospital strains that are treatment-fast, there should be observed current approved medical procedures and techniques. Repeated sterilization or repeated use may constitute a risk for the patient. The tube is exclusively for single use.

Complications

Observance of the current approved medical techniques and performing of a procedure by qualified medical personnel reduces down to a minimum possibility of complications e.g.: Inflammations such as tracheobronchitis, epiglottitis, paresis of hypoglossal or glossal nerves, perforation of swallow, perforation of larynx, pneumothorax, dyspnoea, injuries of lips, tongue, swallow, nose, larynx, vocal apparatus of the larynx, palate, tonsils etc., ulceration of lips, swallow, vocal folds, laryngeal and tracheal stenosis injuries of teeth.

Remarks

The tubes should be stored at a temperature between 5°C and 40°C. Do not use in case of damage of the package. Do not use after the expiry date.

Disposal

After use, dispose of this device in accordance with applicable national regulations on medical waste disposal.

Responsibility – legal guidelines

The Bronchial tubes manufactured by SUMI are subjected to physical and chemical testing and bio-compatibility testing. This allows their safe use and ensures highest quality. Carefully read the instructions for use before using the bronchial tubes. Bronchial tubes made by SUMI are intended for use by qualified medical personnel. The SUMI company will not be liable for the consequences caused by improper selection of tube size or by the use of tubes contrary to the intended purpose and instructions for use.

In the event of a complaint or an occurrence of a medical incident, keep the medical device with its immediate packaging. Every serious medical incident have to be reported to the Manufacturer and National Competent Authority

Size	24Ch	26Ch	28Ch	32Ch	35Ch	37Ch	39Ch	41Ch
Max bronchoscope diameter	2,5	3,4	3,4	4,0	4,0	4,5	4,5	5,0
Suction catheter	6Ch	8Ch	8Ch	10Ch	10Ch	10Ch	10Ch	12Ch

Description

- The bronchial tube is available in performance as the "left-sided" and "right-sided" one.
- It consists of marked bronchial and tracheal tubes and two cuffs, each with a one-way valve and attached pilot balloon.
- The bronchial cuff and corresponding pilot balloon (marked as "bronchial") are in blue colour, whereas the tracheal tube and its balloon (marked as "tracheal") are transparent ones. Connecting tubes are also in blue code colour (bronchial), or transparent (tracheal).
- Printed marks of depth indicate a distance from more distal end of the bronchial segment.
- Each tube has a radiation line for confirmation of correct position of the tube by means of Roentgen rays. Additionally, the "left-sided" version has reference radiographic marks under the cuffs of the form of strips of approx. 20 mm length.
- In the "right sided" hole location in the tube is determined by radiographic line.
- In the version with hook, a hook is additionally fixed on the carina, in the place in which the tube is bent towards bronchus, enabling placing the tube exactly on trachea division.

Indications

Ventilation of a patient during thoracoscopic procedures. Protracted separate bronchial ventilation.

Contraindications

Tracheal stenosis, abnormal anatomy of the airway, sensitization to components of the medical device, the version with a hook to be fixed on the spur may not be used during the surgery performed in the carina area.

Users

Healthcare professionals.

Target group

Patients requiring mechanical ventilation

Instruction for use

- There should be selected appropriate size of the tube individually for each patient.
- To gently remove the sterile bronchial tube from the packaging.
- To firmly mount the 15 mm connectors in the tracheal tube and bronchial junctions so as to prevent them from disconnecting in the middle of using.
- Before use close the bronchoscope port opening on the double swivel elbow connector of the bronchial tube.
- In case of any deformation of the bronchoscope port opening, the special attention should be paid to the leak tightness of the entire system.
- To check the cuff and pilot balloon as well as a valve of each tube by filling before use. To place the end-part of Luer syringe in housing of the inflation valve of the cuff and give an injection of air in the quantity sufficient to completely fill the cuff.
- After pilot inflation of the cuff, deflate it completely.
- To remove a guide and intubate the patient through the mouth or nose paying attention to observance of current approved medical techniques.
- The marker of intubation depth may be of help to you in making corrections of the tube's arrangement.
- To inflate the cuffs, remove the syringe from the valve, quickly connect the tube with a source of ventilation and check correctness of the tube's arrangement, paying attention to observance of current approved medical techniques.
- Check, whether the insufflation system is tight one. Integrity of the inflation system should be monitored both to begin with and periodically during intubation. Pressure in the cuffs should be under precise monitoring and any deviations from the selected pressure should be corrected.
- In order to deflate air of one lung, one should open given lumen of the tube to the room atmosphere prior to mounting a clip on given connecting tube.
- Prior to extubation, each cuff should be deflated through insertion of a syringe in housing of the valve and mixture of gases should be removed through deflation of the pilot balloon.
- Extubate the patient upon observance of the current approved medical techniques.
- Throw the bronchial tube and all accessories aside.

Precautions

- A physician should take appropriate precautions to reduce potential damage of the tracheal cuff as a result of its contact with the patient's teeth or all tools having sharp edges, which constitute a threat for maintaining integrity of the cuff. Therefore, one should proceed with care to avoid damaging of thin walls of the cuff during intubation, which would cause necessity of the patient's extubation and repeated intubation. One should not use a tube having damaged cuff.
- Use of lidocaine spray is associated with formation of pinholes in the PVC cuff.
- Diffusion of the mixture of oxygen and nitrogen, oxygen or air may change volume of the cuff as well as the pressure.
- The cuff should not be over inflated. Ordinarily, the pressure in the cuff should not exceed 25 cm H₂O. Over inflation can result in tracheal damage, rupture of the cuff with subsequent deflation or cuff's distortion, which consequently may bring about airway blockage.
- Overinflating cuff may lead to tracheal damage, deformation or even cuff rupture, what result in cuff deflation or oxygen flow blockage.
- Please be particularly careful in case of introducing the tube with a hook to be fixed on carina.
- Deflate the cuff before rearrangement of the tube. Moving of the tube upon inflated cuff may cause wounding of the patient or damaging of the cuff. After terminated deflation of the cuff, the pilot balloon will be completely collapsed. Check the tube's correct position after each rearrangement thereof.
- Use of the bronchial tube in the procedures, that require using of a laser or other electrosurgical devices with an active electrode in the same place is not allowed.
- When the patient's position or arrangement of the tube undergo a change after intubation, one should check whether the tube is in the correct position. Each change of the tube's position should be immediately corrected.

SUMI sp. z o.o. sp. k.
Drobiarska 35, 05-070 Sulejówek, Poland
tel/fax (+48 22) 783 30 23 www.sumi.com.pl

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