

**NEW**

# MRI-compatible reinforced tracheostomy tube

## SuMRI



**For use in Magnetic Resonance Imaging**



**MEDICAL EQUIPMENT MANUFACTURER**

# SuMRI | MRI-compatible reinforced tracheostomy tube

## Advantages:

Suitability for MRI use.

Insertion **depth marker**.

Innovative reinforcement with a medical-grade plastic spiral.

Made of thermoplastic, **siliconized PVC**:


- **easy insertion**
- **flexibility once placed in the airway**

Control balloon marked with tube size.

Adjustable, soft and flexible flange for regulating tube length.

Tube length adapted for **obese** patients and those with tracheal stenosis.

Soft, **high-volume** low-pressure cylindrical cuff.



## Tube sizes:

Size I.D. [mm]	O.D. [mm]	Tube length range [mm]	Reference number	
			with cuff	without cuff
6,0	8,2	57 - 72	37-6060	38-6060
7,0	9,7	64 - 85	37-7060	38-7060
7,5	10,3	68 - 93	37-7560	38-7560
8,0	11,0	72 - 98	37-8060	38-8060
8,5	11,7	74 - 102	37-8560	38-8560
9,0	12,3	76 - 105	37-9060	38-9060
10,0	13,7	80 - 111	37-1060	38-1060



www.sumi.com.pl



SUMI  
Ul. Drobiarska 35  
05-070 Sulejów



sumi@sumi.com.pl  
(+48 22) 783 30 23