




### Tracheostomy pad




### Heat and moisture exchanger - "artificial nose"



### Tracheostomy tube holder



## Tracheostomy tubes



**SUMINI** tracheostomy tubes were designed with the **safety** and **comfort** of the youngest patients in mind. Made from soft, flexible **medical-grade material**, they are kink-resistant and adapt perfectly to the delicate anatomy of a child. The specially designed flange with internal recesses allows free airflow to the skin of the neck, helping to prevent irritation and pressure sores. The tubes are available in **neonate** and **paediatric** versions, with or without a cuff, and their radiopacity facilitates easy verification of correct positioning under **X-ray**. This combination of innovative technology and attentive care ensures safe and comfortable ventilation in everyday paediatric care.



At the request of parents and physicians, we prepare custom-made tubes tailored individually to the needs of young patients.

**SUMINI** paediatric tracheostomy tube with cuff \*



size I.D. [mm]	O.D. [mm]	length L [mm]	Ø cuff [mm]	ref. number
2,5	4,0	32	8,0	2P-2500
3,0	4,5	35	9,0	2P-3000
3,5	5,1	39	10,0	2P-3500
4,0	5,8	43	12,0	2P-4000
4,5	6,5	47	12,0	2P-4500
5,0	7,1	49	15,0	2P-5000
5,5	7,6	54	15,0	2P-5500
6,0	8,3	60	18,0	2P-6000
6,5	8,9	66	18,0	2P-6500

**SUMINI** neonate tracheostomy tube with cuff \*



size I.D. [mm]	O.D. [mm]	length L [mm]	Ø cuff [mm]	ref. number
2,5	4,0	30	8,0	2N-2500
3,0	4,5	31	9,0	2N-3000
3,5	5,1	33	10,0	2N-3500
4,0	5,8	35	12,0	2N-4000

**SUMINI** paediatric tracheostomy tube without cuff



size I.D. [mm]	O.D. [mm]	length L [mm]	ref. number
2,5	4,0	32	3P-2500
3,0	4,5	35	3P-3000
3,5	5,1	39	3P-3500
4,0	5,8	43	3P-4000
4,5	6,5	47	3P-4500
5,0	7,1	49	3P-5000
5,5	7,6	54	3P-5500
6,0	8,3	60	3P-6000
6,5	8,9	66	3P-6500

**SUMINI** neonate tracheostomy tube without cuff



size I.D. [mm]	O.D. [mm]	length L [mm]	ref. number
2,5	4,0	30	3N-2500
3,0	4,5	31	3N-3000
3,5	5,1	33	3N-3500
4,0	5,8	35	3N-4000

Advantages of **SUMINI** tubes:

- made of **high-quality**, flexible medical material that adapts to patient anatomy,
- high-volume **low-pressure** cuff,\*
- soft and **flexible flange** that easily conforms to the shape of the neck,
- specially designed flange ensures secure fixation and minimizes the risk of skin irritation,
- perfectly **matched obturator** facilitating easy tube insertion,
- neck **tape** for fastening,
- visible in **X-rays**,
- rounded, **atraumatic tip**,
- latex-free, **phthalate-free**,
- **sterile**, single-use.



To ensure the best possible fit to the patient's anatomy, depending on the type of flange, the tubes have different lengths for the same size. Recommended for long-term ventilation.